

A CALIBRE SCIENTIFIC COMPANY





Reagecon Buffer Standard pH 12.00 at 25°C ISO 17034 Certified Reference Material

Item No. CRM112025

Reagecon's Buffer Standard pH 12.00 at 25°C Certified Reference Material is produced under our ISO 17034 (A2LA Ref: 6739.01) Accreditation and is NIST traceable.

Additional Details

Pack Size 1L
Weight 1.16
Weight Unit Kg

Tariff Code 3822 190090

Country of

Origin Ireland

DESCRIPTION

TRACEABILITY: Directly traceable to the IUPAC pH scale by an unbroken chain of traceability. Reagecon achieve this
traceability through a series of comparisons with the key reference materials being Standard Reference Materials (SRMs)
manufactured by NIST. Reagecon's pH Buffer Standards meet the ISO definition of metrological traceability "The ability to
relate measurements back to a stated reference (usually an international standard) through an unbroken chain of
comparisons each having stated uncertainties of measurement." Reagecon's traceability claims are independently verified
through our accreditation to ISO 17034 and ISO 17025.

- ISO 17034 accreditation requirements include raw material characterisation; homogeneity and stability studies; assignmen of certified values as well as packaging and documentation control.
- STABILITY AND SHELF LIFE: Manufactured to exacting specifications with an extended shelf life and stability; even after
 opening the bottle thereby eliminating the requirement to open a fresh bottle of Standard every time it is used.
- TEMPERATURE DEPENDENCE & EXPIRY: Temperature dependence data is available via the product label as are lot numbers and expiry dates for user convenience.
- CALIBRATION BUFFERS: Reagecon pH Buffers are pre-programmed into the instruments of most major manufacturers.
- CONTROL BUFFERS: For increased confidence in their test measurements; analysts should regularly measure the pH of a
 Control Standard. If an acceptable value is obtained from the Control Standard measurement then the analysts can have
 improved confidence that their test measurements will be correct. Reagecon's extensive range of pH Buffers means that
 there will be a Reagecon pH Buffer which can be used as a Control Buffer for all pH applications.
- WHY USE TRACEABLE pH BUFFERS: Your pH measurements can only be as good as the pH Buffers that you use. If your pH
 calibration is made using traceable pH Buffers then you have a direct link to the International pH scale for your
 measurements. Without this link you are not entitled to report your measurements in pH units so the number displayed on
 your pH meter is just that a number and is not a pH value. Traceable analysis is necessary for consistency and universal
 acceptance of your pH results including acceptance by regulatory bodies.

PRODUCT SPECIFICATIONS

pH Buffer Reference Standard; pH Buffer Reference Material

Value	pH 12.00 ± 0.05
Additional Information	Specification ± 0.05 pH units
Tested	at 25°C

Use	Laboratory Chemical
Packaging	HDPE Bottle
Physical Form	Liquid